



DOCTOR'S PREFERENCE GUIDE

Please complete the questionnaire regarding general preferences and designs.
Thank you!

DOCTOR _____
 CLINIC _____
 ADDRESS _____
 CITY/ZIP _____
 BUSINESS PHONE _____
 BUSINESS FAX _____
 HOME PHONE _____
 E-MAIL _____
 ASSISTANTS _____
 HYGIENISTS _____
 RECEPTIONISTS _____

OFFICE HOURS

MON	TUE	WED	THUR	FRI	SAT

LUNCH HOURS

MON	TUE	WED	THUR	FRI	SAT

OFFICE CLOSED DURING LUNCH YES _____ NO _____

DROP SITE YES _____ NO _____

LOCATION _____

MAILING STICKERS YES _____ NO _____

DIE TRIMMING ■ ■ ■ ■ ■ ■

- Trim to visible margin.
- Trim past visible margin, until tooth structure is no longer apparent or an undercut occurs.
- Trim past visible margin to extend margin 1 mm.
- Return to doctor for die trim.

DIE SPACER ■ ■ ■ **RIDGE RELIEF** ■ ■ ■

- | | |
|----------------------------------|---------------------------------|
| <input type="checkbox"/> None | <input type="checkbox"/> None |
| <input type="checkbox"/> 1 Coat | <input type="checkbox"/> Slight |
| <input type="checkbox"/> 2 Coats | <input type="checkbox"/> Medium |
| <input type="checkbox"/> 3 Coats | <input type="checkbox"/> Heavy |
| <input type="checkbox"/> 4 Coats | |

PORCELAIN ALLOYS ■ ■ ■ ■ ■ ■

- High Gold Yellow (High Noble)
- No Gold White (Noble)
- High Gold White (High Noble)

FULL CAST ALLOYS ■ ■ ■ ■ ■ ■

- High Gold Yellow (High Noble)
- High Gold White (High Noble)
- No Gold White (Noble)

PORCELAIN DESIGN ■ ■ ■ ■ ■ ■

Check your preference:

- Metal linguals instead of an opaque appearance
- Occlusal islands instead of opaque anatomy
- Opaque linguals instead of bulky linguals
- Have lab call for opposing tooth reduction


Comments:

ANTERIORES _____

■ ■ ■ ■ ■ ■

Circle the designs you prefer:

PORCELAIN DESIGN



- Porcelain Shoulder
- 360° Metal Margin
- Metal Contacts
- Metal Occlusal Island

Comments:




POSTERIORES _____

PONTIC DESIGN ■ ■ ■ ■ ■



Circle the design you prefer:



ANTERIOR INTERPROXIMALS ■ ■ ■



POSTERIOR INTERPROXIMALS ■ ■ ■



COMMENTS ON LINGUAL ANATOMY AND DEVELOPMENT OF CINGULUM AS WELL AS BUCCAL-LINGUAL WIDTH OF PONTIC:

MESIAL AND DISTAL CONTACT ■ ■ ■


- Match dental arch
- Light broad contact
- Firm broad contact

OCCLUSAL CONTACTS ■ ■ ■ ■ ■

- Foil
- No foil
- Space open in centric

OCCLUSAL STAIN ■ ■ ■ ■ ■

- None Orange
- Brown Other
- Match shade tab (follow tendencies)



GLAZE ■ ■ ■ ■ ■

- Low (dull satin finish)
- Medium (smooth shiny surface) as shade tab
- High (glossy finish)

BRAND OF SHADE TAB ■ ■ ■ ■ ■

BRAND OF CEMENT ■ ■ ■ ■ ■

BRAND OF IMPRESSION MATERIAL ■ ■ ■ ■ ■

COMMENTS ■ ■ ■ ■ ■

Thank you for your time and cooperation.